

BONE CURETTES



SECTION H



To scrape, shape and clean bone.

CONTENTS

Bone Curettes	1-6
Endoscopic Bone Curettes	7-8
Long Bone Curettes	9-10
Micro Curettes	11-12



OM 29-0632
6" (152mm)
double ended
oval cups

JANSEN CURETTE




OM 29-0620
5 1/2" (140mm)
double ended
4mm and 5mm round cups

MARTINI BONE CURETTE



OM 29-0630
5 1/2" (140mm)
double ended
3mm and 4mm oval cups


WILLIGER BONE CURETTE



OM 46-0560 5" double ended
OM 46-0562 6 1/2" double ended
OM 46-0561 8" double ended

VOLKMANN CURETTE
Oval and round cups

BONE AWLS




OM 49-0200 OM 49-0202

OM 46-0550 OM 49-0204

OM 49-0200 5" (127mm)
 OM 49-0202 5 1/2" (140mm)
 OM 46-0550 6 1/2" (165mm)
 OM 49-0204 8" (203mm)

VOLKMANN CURETTE

Oval cups, double ended



3mm
4mm
5mm

3mm
4mm
5mm

3mm
4mm
5mm

3mm
4mm
5mm

Vertical angle 45°
 OM 49-0403 3mm diameter
 OM 49-0404 4mm diameter
 OM 49-0405 5mm diameter

Vertical angle 90°
 OM 49-0413 3mm diameter
 OM 49-0414 4mm diameter
 OM 49-0415 5mm diameter

Horizontal angle 45°
 OM 49-0423 3mm diameter
 OM 49-0424 4mm diameter
 OM 49-0425 5mm diameter

Horizontal angle 90°
 OM 49-0433 3mm diameter
 OM 49-0434 4mm diameter
 OM 49-0435 5mm diameter

RAY PITUITARY CURETTES

7 1/2"





OM 49-0580 8" (203mm) #000
 OM 49-0582 8" (203mm) #00
 OM 49-0584 8" (203mm) #0

EPSTEIN CURETTE
 With hollow handle



OM 49-0574 straight
 OM 49-0576 down angle
 OM 49-0578 up angle

SCOVILLE CURETTE
 With hollow handle
 10" (254mm)




6 3/4" (172mm) straight
 OM 49-0120 #3/0
 OM 49-0121 #2/0
 OM 49-0122 #0
 OM 49-0123 #1
 OM 49-0124 #2
 OM 49-0125 #3
 OM 49-0126 #4



6 3/4" (172mm) angled 45°
 OM 49-0120A #3/0
 OM 49-0121A #2/0
 OM 49-0122A #0
 OM 49-0123A #1
 OM 49-0124A #2
 OM 49-0125A #3
 OM 49-0126A #4

BRUN (SPRATT) CURETTE
 With hollow handle, round cup

BONE CURETTES



7" (178mm) straight
 OM 46-0511 #5/0
 OM 46-0512 #4/0
 OM 46-0513 #3/0
 OM 46-0515 #2/0
 OM 46-0517 #0
 OM 46-0519 #1
 OM 46-0521 #2
 OM 46-0523 #3
 OM 46-0525 #4
 OM 46-0529 #5
 OM 46-0531 #6


7" (178mm) angled 45°
 OM 46-0511A #5/0
 OM 46-0512A #4/0
 OM 46-0513A #3/0
 OM 46-0515A #2/0
 OM 46-0517A #0
 OM 46-0519A #1
 OM 46-0521A #2
 OM 46-0523A #3
 OM 46-0525A #4
 OM 46-0529A #5
 OM 46-0531A #6

9" (229mm) straight
 OM 46-0563 #3/0
 OM 46-0564 #2/0
 OM 46-0565 #0
 OM 46-0566 #1
 OM 46-0567 #2
 OM 46-0568 #3
 OM 46-0569 #4
 OM 46-0570 #5
 OM 46-0571 #6

9" (229mm) angled 45°
 OM 49-0009 #4/0
 OM 49-0110 #3/0
 OM 49-0111 #2/0
 OM 49-0112 #0
 OM 49-0113 #1
 OM 49-0114 #2
 OM 49-0115 #3
 OM 49-0116 #4
 OM 49-0117 #5
 OM 49-0118 #6

BRUN (SPRATT) CURETTE

With hollow handle, oval cup



Straight
 OM 46-0551 #3/0
 OM 46-0552 #2/0
 OM 46-0553 #0
 OM 46-0554 #1
 OM 46-0555 #2
 OM 46-0556 #3
 OM 46-0557 #4
 OM 46-0558 #5
 OM 46-0559 #6

Angled 45°
 OM 49-0150 #3/0
 OM 49-0151 #2/0
 OM 49-0152 #0
 OM 49-0153 #1
 OM 49-0154 #2
 OM 49-0155 #3
 OM 49-0156 #4
 OM 49-0157 #5
 OM 49-0158 #6

BRUN CURETTE

9" (229mm) with solid hexagonal handle



Anodized aluminum handle
stainless steel, working end

All stainless.

Straight

- OM 46-0422 #3/0
- OM 46-0423 #2/0
- OM 46-0424 #0
- OM 46-0425 #1
- OM 46-0426 #2
- OM 46-0427 #3
- OM 46-0428 #4
- OM 46-0429 #5
- OM 46-0430 #6

Angled

- OM 46-0433 #3/0
- OM 46-0434 #2/0
- OM 46-0435 #0
- OM 46-0436 #1
- OM 46-0437 #2
- OM 46-0438 #3
- OM 46-0439 #4
- OM 46-0440 #5
- OM 46-0441 #6

Straight

- OM 46-0422S #3/0
- OM 46-0423S #2/0
- OM 46-0424S #0
- OM 46-0425S #1
- OM 46-0426S #2
- OM 46-0427S #3
- OM 46-0428S #4
- OM 46-0429S #5
- OM 46-0430S #6

Angled

- OM 46-0433S #3/0
- OM 46-0434S #2/0
- OM 46-0435S #0
- OM 46-0436S #1
- OM 46-0437S #2
- OM 46-0438S #3
- OM 46-0439S #4
- OM 46-0440S #5
- OM 46-0441S #6

COBB CURETTE

11" (280mm) with knurled handle



- OM PS03.959-0 Cobb Curette 11 1/4" Straight #0
- OM PS03.959-00 Cobb Curette 11 1/4" Straight #00
- OM PS03.959-000 Cobb Curette 11 1/4" Straight #000
- OM PS03.959-1 Cobb Curette 11 1/4" Straight #1
- OM PS03.959-2 Cobb Curette 11 1/4" Straight #2
- OM PS03.959-3 Cobb Curette 11 1/4" Straight #3
- OM PS03.959-4 Cobb Curette 11 1/4" Straight #4
- OM PS03.959-5 Cobb Curette 11 1/4" Straight #5
- OM PS03.959-6 Cobb Curette 11 1/4" Straight #6

SUNY UPSTATE CURETTE

11 1/4" (286mm)



Silicone handles are available in assorted colors per request.


OM W45-234
Sharp, serrated with black silicone handle left

OM W45-235
Sharp, serrated with black silicone handle right

OM W45-233
Sharp, serrated with black silicone handle straight

CURETTE #8

BONE CURETTES



Straight

OM 49-0450 size #1, 8mm diameter
 OM 49-0452 size #2, 6mm diameter
 OM 49-0454 size #3, 3mm diameter

Angled 35°

OM 49-0456 size #1, 8mm diameter
 OM 49-0458 size #2, 6mm diameter
 OM 49-0460 size #3, 3mm diameter

CONE RING CURETTE

9" (229mm) with knurled handle with fenestrated cup



OM 55-0680 10" (250mm)
 OM 55-0684 11 1/2" (292mm)



OM 55-0681 10" (250mm)
 OM 55-0685 11 1/2" (292mm)



OM 55-0682 10 1/4" (260mm)
 OM 55-0686 12" (305mm)

RING CURETTE

Hollow handle with fenestrated cup




OM 29-0503



OM 29-0505



OM 29-0507



OM 29-0495 1mm
 OM 29-0496 2mm
 OM 29-0497 3mm
 OM 29-0498 4mm
 OM 29-0499 5mm
 OM 29-0500 6mm



PIFFARD CURETTE

5 1/2" (140mm)




FOX CURETTE


5 1/2" (140mm) round




OM 55-0532 17" (432mm)
40° angle
2mm x 4mm cup
2.7" tip length
9.3" shaft length
5" handle length
.285" (7.2mm) shaft diameter




OM 55-0523 17" (432mm)
straight
3mm x 6mm cup
2.7" tip length
9.3" shaft length
5" handle length
.285" (7.2mm) shaft diameter



OM 55-0533 17" (432mm)
40° angle
3mm x 6mm cup
2.7" tip length
9.3" shaft length
5" handle length
.285" (7.2mm) shaft diameter





OM 55-0543 20" (508mm)
40° angle
3mm x 6mm cup
3.5" tip length
10" shaft length
6.5" handle length
.45" (11.5mm) shaft diameter




OM 55-0556 20" (508mm)
30° angle
6mm x 9mm cup
3.5" tip length
10" shaft length
6.5" handle length
.45" (11.5mm) shaft diameter

Phenolic handles are autoclavable to 290°F (143°C).

ENDOSCOPIC BRUN CURETTE
Graduated with phenolic handle




OM 55-0593 17" (432mm)
reverse 40° angle
3mm x 6mm cup
2.7" tip length
9.3" shaft length
5" handle length
.285" (7.2mm) shaft diameter



OM 55-0586 20" (508mm)
reverse 30° angle
6mm x 9mm cup
3.5" tip length
10" shaft length
6.5" handle length
.45" (11.5mm) shaft diameter

Phenolic handles are autoclavable to 290°F (143°C).



ENDOSCOPIC EPSTEIN CURETTE
Graduated with phenolic handle




OM 55-0510 20" (508mm)
straight
9mm x 15mm cup
3.5" tip length
10" shaft length
6.5" handle length
.45" (11.5mm) shaft diameter

Phenolic handles are autoclavable to 290°F (143°C).

ENDOSCOPIC SCOVILLE CURETTE
With phenolic handle



OM 55-0573 17" (432mm)
40° angle
3mm cup diameter
2.7" tip length
9.3" shaft length
5" handle length
.285" (7.2mm) shaft diameter



OM 55-0566 20" (508mm)
30° angle
6mm cup diameter
3.5" tip length
10" shaft length
6.5" handle length
.45" (11.5mm) shaft diameter

Phenolic handles are autoclavable to 290°F (143°C).

ENDOSCOPIC CONE RING CURETTE
With phenolic handle



Straight 15" (381mm)
 OM PS00 991-01 #4/0
 OM PS00 991-02 #3/0
 OM PS00 991-03 #2/0
 OM PS00 991-04 #0
 OM PS00 991-05 #1
 OM PS00 991-06 #2
 OM PS00 991-07 #3
 OM PS00 991-08 #4
 OM PS00 991-09 #5
 OM PS00 991-10 #6

Angled 15" (381mm)
 OM PS00 991-01A #4/0
 OM PS00 991-02A #3/0
 OM PS00 991-03A #2/0
 OM PS00 991-04A #0
 OM PS00 991-05A #1
 OM PS00 991-06A #2
 OM PS00 991-07A #3
 OM PS00 991-08A #4
 OM PS00 991-09A #5
 OM PS00 991-10A #6

LONG BONE CURETTE



Straight 15" (381mm)
 OM 55-0576 (10mm X 19mm cup)
 OM 55-0577 (13mm X 25mm cup)
 OM 55-0578 (15mm X 29mm cup)

VOLKMANN CURETTE

Straight

OM 93-0112 15" (381mm) overall, #3/0
OM 93-0113 15" (381mm) overall, #2/0
OM 93-0114 15" (381mm) overall, #0
OM 93-0115 15" (381mm) overall, #1
OM 93-0116 15" (381mm) overall, #2
OM 93-0117 15" (381mm) overall, #3
OM 93-0118 15" (381mm) overall, #4
OM 93-0119 15" (381mm) overall, #5
OM 93-0125 15" (381mm) overall, #6
OM 93-0127 15" (381mm) overall, #7



LONG COBB CURETTE
15" overall length with 9" handle

Angled

OM 93-0205 15" (381mm) overall, #3/0
OM 93-0206 15" (381mm) overall, #2/0
OM 93-0207 15" (381mm) overall, #0
OM 93-0208 15" (381mm) overall, #1
OM 93-0209 15" (381mm) overall, #2
OM 93-0210 15" (381mm) overall, #3
OM 93-0211 15" (381mm) overall, #4
OM 93-0212 15" (381mm) overall, #5
OM 93-0213 15" (381mm) overall, #6
OM 93-0214 15" (381mm) overall, #7



OM RD 11-1142 size 6
OM RD 11-1143 size 8
OM RD 11-1144 size 10
OM RD 11-1145 size 12
OM RD 11-1146 size 15
OM RD 11-1147 size 18

MCELROY CURETTE
14" (355mm) stainless steel

CERVICAL MICRO CURETTE

Forward Straight

OM 99-0001 #0
 OM 99-0002 #00
 OM 99-0003 #000
 OM 99-0004 #0000
 OM 99-0005 Bone Cutter #00

Forward Angled

OM 99-0006 #0
 OM 99-0007 #00
 OM 99-0008 #000
 OM 99-0009 #0000
 OM 99-0010 Bone Cutter #00

Backward Straight

OM 99-0011 #0
 OM 99-0012 #00
 OM 99-0013 #000
 OM 99-0014 #0000
 OM 99-0015 Bone Cutter #00

Backward Angle

OM 99-0016 #0
 OM 99-0017 #00
 OM 99-0018 #000
 OM 99-0019 #0000
 OM 99-0020 Bone Cutter #00

OM 99-0021 Cervical Micro Nerve Hook
 OM 99-0022 Cervical Ultra-Micro Nerve Hook
 OM 99-0023 Cervical Micro Waxer
 OM 99-0024 Cervical Micro Knife Handle

XL CERVICAL MICRO CURETTE

Forward Straight

OM 99-0030 #0
 OM 99-0031 #00
 OM 99-0032 #000
 OM 99-0033 #0000
 OM 99-0034 Bone Cutter

Forward Angled

OM 99-0035 #0
 OM 99-0036 #00
 OM 99-0037 #000
 OM 99-0038 #0000
 OM 99-0039 Bone Cutter

Backward Straight

OM 99-0040 #0
 OM 99-0041 #00
 OM 99-0042 #000
 OM 99-0043 #0000
 OM 99-0044 Bone Cutter

Backward Angle

OM 99-0045 #0
 OM 99-0046 #00
 OM 99-0047 #000
 OM 99-0048 #0000
 OM 99-0049 Bone Cutter

OM 99-0050 XL Cervical Micro Nerve Hook
 OM 99-0051 XL Cervical Ultra-Micro Nerve Hook
 OM 99-0052 Cervical Micro Waxer
 OM 99-0053 Cervical Micro Knife Handle



LUMBAR MICRO CURETTE

Forward Straight

OM 99-0060 #0
 OM 99-0061 #1
 OM 99-0062 #2

Forward Angled

OM 99-0063 #0
 OM 99-0064 #1
 OM 99-0065 #2

Forward Down

OM 99-0066 #0
 OM 99-0067 #1
 OM 99-0068 #2

Backward Straight

OM 99-0069 #0
 OM 99-0070 #1
 OM 99-0071 #2

MICRO CURETTES

OM 99-0300	Sterilizing Tray with Silicone Holders for micro curettes		
OM 99-0301	Microsect curette	25.9cm Fig. 1	forward angle
OM 99-0302	Microsect curette	25.9cm Fig. 1	backward angle
OM 99-0303	Microsect curette	25.9cm Fig. 2	forward angle
OM 99-0304	Microsect curette	25.9cm Fig. 2	backward angle
OM 99-0305	Microsect curette	25.9cm Fig. 3	forward angle
OM 99-0306	Microsect curette	25.9cm Fig. 3	backward angle
OM 99-0307	Microsect curette	25.9cm Fig. 4	forward angle
OM 99-0308	Microsect curette	25.9cm Fig. 4	backward angle
OM 99-0309	Microsect curette	25.9cm Fig. 5	forward angle
OM 99-0310	Microsect curette	25.9cm Fig. 5	backward angle
OM 99-0311	Bray End Dissector	29.5cm	5mm
OM 99-0312	Bray End Dissector	29.5cm	3mm



MICRO CURETTES