



# OSTEOTOMES, CHISELS AND GOUGES



## SECTION J



To score, cut, scrape, clean and sculpt bone.

### CONTENTS

Interchangeable Set	1
Chisels	4-9
Gouges	5-12
Osteotomes	2-12





- OM 49-0668 5mm
- OM 49-0672 10mm
- OM 49-0676 16mm
- OM 49-0680 25mm

**OSTEOTOME BLADES**




OM 49-0672



OM 49-0664

- OM 49-0652 5mm
- OM 49-0656 10mm
- OM 49-0660 16mm
- OM 49-0664 25mm

**CHISEL BLADES**




**KEY**  
OM 49-0634

Phenolic handle is autoclavable to 290°F (143°C).

- OM 49-0630 complete set with handle, key and all blades
- OM 49-0632 7" black phenolic handle handpiece only

**INTERCHANGEABLE OSTEOTOME, CHISEL AND GOUGE SET**



OM 49-0640

OM 49-0644


OM 49-0648


- OM 49-0640 1" width, 50mm radius
- OM 49-0644 1" width, 60mm radius
- OM 49-0648 1" width, 70mm radius


Blades should be discarded after use.


**GOUGE BLADES**


# OSTEOTOMES, CHISELS AND GOUGES





  
4mm

  
6mm

  
8mm


  
10mm


  
12mm


  
15mm


5" Straight			
OM 46-0337	2mm	OM 46-0343	14mm
OM 46-0336	4mm	OM 46-0345	15mm
OM 46-0338	6mm	OM 46-0346	18mm
OM 46-0340	8mm	OM 46-0347	19mm
OM 46-0342	10mm	OM 46-0335	20mm
OM 46-0344	12mm		

**MINI-LAMBOTTE OSTEOTOME**





  
10mm

  
20mm

The unusually long bevel allows the cut bone to bend gradually as it is cut.

7" Straight With calibration lines		7" Curved With calibration lines	
OM 69-0940	10mm	OM 69-0970	5mm
OM 69-0946	20mm	OM 69-0975	10mm
OM 69-0950	25mm		

**LONG BEVEL OSTEOTOME**



  
4mm

  
8mm

  
10mm


  
12mm

  
25mm

7" Straight With calibration lines		7" Curved With calibration lines	
OM 46-0348	4mm	OM 46-0348C	4mm
OM 46-0349	5mm	OM 46-0349C	5mm
OM 46-0350	6mm	OM 46-0350C	6mm
OM 46-0358	8mm	OM 46-0358C	8mm
OM 46-0351	10mm	OM 46-0351C	10mm
OM 46-0352	12mm	OM 46-0352C	12mm
OM 46-0353	16mm	OM 46-0353C	16mm
OM 46-0354	19mm	OM 46-0354C	19mm
OM 46-0355	22mm	OM 46-0355C	22mm
OM 46-0356	25mm	OM 46-0356C	25mm

**LAMBOTTE OSTEOTOME**

Racks for 9" Lambotte  
Osteotomes are  
shown on page AA2.



**9" Straight**  
Without calibration lines


OM 46-0304	1/4"	6mm
OM 46-0306	1/2"	13mm
OM 46-0310	3/4"	19mm
OM 46-0314	1"	25mm
OM 46-0318	1 1/4"	32mm
OM 46-0322	1 1/2"	38mm

**9" Curved**  
Without calibration lines

OM 46-0328	1/4"	6mm
OM 46-0329	1/2"	13mm
OM 46-0330	3/4"	19mm
OM 46-0331	1"	25mm
OM 46-0332	1 1/4"	32mm
OM 46-0333	1 1/2"	38mm

## LAMBOTTE OSTEOTOME

Racks for Hoke Osteotomes  
are shown on page AA3.



**Straight**

OM 46-0497	1/8"	3mm
OM 46-0498	3/16"	4mm
OM 46-0499	1/4"	6mm
OM 46-0500	5/16"	8mm
OM 46-0501	3/8"	10mm
OM 46-0502	1/2"	13mm
OM 46-0503	5/8"	16mm
OM 46-0504	3/4"	19mm

**Curved**


OM 46-0700	1/8"	3mm
OM 46-0701	3/16"	4mm
OM 46-0702	1/4"	6mm
OM 46-0703	5/16"	8mm
OM 46-0704	3/8"	10mm
OM 46-0706	1/2"	13mm
OM 46-0708	5/8"	16mm
OM 46-0709	3/4"	19mm

## HOKE OSTEOTOME

5 1/4" (140mm) straight


# OSTEOTOMES, CHISELS AND GOUGES

Racks for Hoke Chisels are shown in Section AA.




## HOKE CHISEL

5 1/4" (133mm)  
straight



## SHEEHAN OSTEOTOME

6 1/4" (159mm)  
straight




OM 58-0102      OM 58-0104

OM 58-0100 straight  
OM 58-0102 curved, guard right  
OM 58-0104 curved, guard left

## SILVER OSTEOTOME

7" (178mm)  
with single guard




OM 46-0284      6mm  
OM 46-0286      12mm  
OM 46-0288      18mm

## U.S. ARMY PATTERN OSTEOTOME

7" (178mm)



6mm



12mm

OM 46-0289 6mm  
 OM 46-0290 12mm  
 OM 46-0291 14mm

**U.S. ARMY PATTERN CHISEL**

7" (178mm)



6mm



12mm



14mm

OM 46-0292 6mm  
 OM 46-0293 12mm  
 OM 46-0294 14mm

**U.S. ARMY PATTERN GOUGE**

7" (178mm)



4mm




6mm




8mm



10mm



12mm



14mm

OM 49-0040 4mm  
 OM 49-0041 6mm  
 OM 49-0042 8mm  
 OM 49-0043 10mm  
 OM 49-0044 12mm  
 OM 49-0045 14mm

**ALEXANDER CHISEL**

7" (178mm)



4mm




6mm




8mm



10mm



12mm




14mm

OM 49-0050 4mm  
 OM 49-0051 6mm  
 OM 49-0052 8mm  
 OM 49-0053 10mm  
 OM 49-0054 12mm  
 OM 49-0055 14mm

**ALEXANDER GOUGE**


7" (178mm)



- OM 49-0260 4mm
- OM 49-0262 6mm
- OM 49-0264 8mm
- OM 49-0266 10mm
- OM 49-0268 12mm

**LEXER MINI CHISEL**

7" (178mm)



- OM 49-0246 5mm
- OM 49-0248 7mm
- OM 49-0250 10mm
- OM 49-0252 15mm
- OM 49-0254 20mm
- OM 49-0256 25mm
- OM 49-0258 30mm

**LEXER CHISEL**

8 3/4" (222mm)

Racks for Smith Peterson Osteotomes are shown in Section AA.



1/4"



1/2"



1"



**8" Straight**

OM 46-0448	1/4"	6mm
OM 46-0449	3/8"	10mm
OM 46-0450	1/2"	13mm
OM 46-0451	5/8"	16mm
OM 46-0452	3/4"	19mm
OM 46-0453	1"	25mm
OM 46-0454	1 1/4"	32mm

**8" Curved**

OM 46-0457	1/4"	6mm
OM 46-0458	3/8"	10mm
OM 46-0459	1/2"	13mm
OM 46-0460	5/8"	16mm
OM 46-0461	3/4"	19mm
OM 46-0462	1"	25mm
OM 46-0463	1 1/4"	32mm

**SMITH PETERSON OSTEOTOME**

Racks for Smith Peterson Gouges are shown in Section AA.



1/4"



1/2"



3/4"



**8" Straight**

OM 46-0466	1/4"	6mm
OM 46-0467	3/8"	10mm
OM 46-0468	1/2"	13mm
OM 46-0476	5/8"	16mm
OM 46-0477	3/4"	19mm
OM 46-0478	1"	25mm
OM 46-0479	1 1/4"	32mm

**8" Curved**

OM 46-0487	1/4"	6mm
OM 46-0488	3/8"	10mm
OM 46-0489	1/2"	13mm
OM 46-0490	5/8"	16mm
OM 46-0491	3/4"	19mm
OM 46-0493	1"	25mm
OM 46-0495	1 1/4"	32mm

**SMITH PETERSON GOUGE**



# OSTEOTOMES, CHISELS AND GOUGES



OM 49-0002 10mm  
 OM 49-0004 14mm  
 OM 49-0006 18mm  
 OM 49-0008 22mm

**STILLE-TYPE OSTEOTOME**  
 8" (203mm) straight



OM 49-0010 10mm  
 OM 49-0012 14mm  
 OM 49-0014 18mm  
 OM 49-0016 22mm

**STILLE-TYPE CHISEL**  
 8" (203mm) straight



OM 49-0018 10mm  
 OM 49-0020 14mm  
 OM 49-0022 18mm  
 OM 49-0024 22mm

**STILLE-TYPE GOUGE**  
 8" (203mm) straight

Racks for Hibbs Osteotomes are shown in Section AA.



1"



1 1/4"



**9" (229mm) Straight**

OM 46-0359 1/4" (6mm)  
 OM 46-0360 3/8" (10mm)  
 OM 46-0361 1/2" (13mm)  
 OM 46-0362 5/8" (16mm)  
 OM 46-0363 3/4" (19mm)  
 OM 46-0364 7/8" (22mm)  
 OM 46-0365 1" (25mm)  
 OM 46-0366 1 1/8" (29mm)  
 OM 46-0367 1 1/4" (31mm)  
 OM 46-0368 1 1/2" (38mm)


**9" (229mm) Curved**

OM 46-0373 1/4" (6mm)  
 OM 46-0376 3/8" (10mm)  
 OM 46-0377 1/2" (13mm)  
 OM 46-0378 5/8" (16mm)  
 OM 46-0379 3/4" (19mm)  
 OM 46-0380 7/8" (22mm)  
 OM 46-0381 1" (25mm)  
 OM 46-0382 1 1/8" (29mm)  
 OM 46-0383 1 1/4" (31mm)  
 OM 46-0384 1 1/2" (38mm)

**HIBBS OSTEOTOME**



Racks for Hibbs Chisel are shown in Section AA.



9" (229mm) Straight

OM 49-0059	1/4" (6mm)
OM 49-0060	3/8" (10mm)
OM 49-0061	1/2" (13mm)
OM 49-0062	5/8" (16mm)
OM 49-0063	3/4" (19mm)
OM 49-0064	7/8" (22mm)
OM 49-0065	1" (25mm)
OM 49-0066	1 1/8" (29mm)
OM 49-0067	1 1/4" (31mm)
OM 49-0068	1 1/2" (38mm)


9" (229mm) Curved

OM 49-0073	1/4" (6mm)
OM 49-0074	3/8" (10mm)
OM 49-0075	1/2" (13mm)
OM 49-0076	5/8" (16mm)
OM 49-0077	3/4" (19mm)
OM 49-0078	7/8" (22mm)
OM 49-0079	1" (25mm)
OM 49-0080	1 1/8" (29mm)
OM 49-0081	1 1/4" (31mm)
OM 49-0082	1 1/2" (38mm)

**HIBBS CHISEL**



Racks for Hibbs Gouge are shown in Section AA.




9" (229mm) Straight

OM 46-0387	1/4" (6mm)
OM 46-0388	3/8" (10mm)
OM 46-0389	1/2" (13mm)
OM 46-0390	5/8" (16mm)
OM 46-0391	3/4" (19mm)
OM 46-0392	7/8" (22mm)
OM 46-0393	1" (25mm)
OM 46-0394	1 1/8" (29mm)
OM 46-0395	1 1/4" (31mm)
OM 46-0401	1 1/2" (38mm)


9" (229mm) Curved

OM 46-0402	1/4" (6mm)
OM 46-0403	3/8" (10mm)
OM 46-0404	1/2" (13mm)
OM 46-0405	5/8" (16mm)
OM 46-0407	3/4" (19mm)
OM 46-0408	7/8" (22mm)
OM 46-0409	1" (25mm)
OM 46-0410	1 1/8" (29mm)
OM 46-0411	1 1/4" (31mm)
OM 46-0412	1 1/2" (38mm)


**HIBBS GOUGE**



**SHOULDER PENETRATING GOUGE**



**SWAN NECK GOUGE**




**COBB OSTEOTOME**



OM 49-0030  
straight  
5/32" (7mm) blade width



OM 49-0032  
medium curve  
5/32" (7mm) blade width



OM 49-0034  
strong curve  
5/32" (7mm) blade width




OM 49-0036  
reverse curve  
5/32" (7mm) blade width

**COBB-TYPE SPINAL GOUGE**


11" (280mm)  
All stainless, with knurled handle

Racks for Cobb-Type Gouges  
are shown in Section AA.



OM 55-0630  
18 1/2" (470mm) with stainless steel handle  
.28" (7mm) tip width straight  
3 1/2" tip length  
10 3/4" shaft length  
4 1/4" handle length  
.45" (11.5mm) shaft diameter

**ENDOSCOPIC HIBBS GOUGE**



OM 55-0642  
18 1/2" (470mm) with stainless steel handle  
3/8" (10mm) tip width straight  
3 1/2" tip length  
10 3/4" shaft length  
4 1/4" handle length  
.45" (11.5mm) shaft diameter

**ENDOSCOPIC HIBBS OSTEOTOME**