

# MEASURING DEVICES

## SECTION W



To determine dimensions during surgery.


### CONTENTS

Screw Depth Gauge	1
Calipers	1-2
Rulers	3
Easy Clean Depth Gauge	4



**TOWNLEY CALIPER**

OM 48-0314  
 4" (102mm)  
 takes inside and outside measurements up to 4" (10cm)



**BONE COMPASS**

OM 48-0302  
 8" (203mm)

*Used to determine screw length required during internal fracture fixation (mm graduated). Gauges snap apart for easy cleaning.*



**MINI SCREW DEPTH GAUGE**

OM 47-0934  
 4 3/4" (121mm)  
 up to 30mm



**SMALL SCREW DEPTH GAUGE**

OM 47-0936  
 6 1/4" (159mm)  
 up to 50mm




**LARGE SCREW DEPTH GAUGE**

OM 47-0940  
 8 3/4" (222mm)  
 up to 100mm




**X-LARGE SCREW DEPTH GAUGE**

OM 47-0935  
 11" (279mm)  
 up to 120mm




**CASTROVIEJO CALIPER**

OM 53-0040  
3 1/2" (89mm)  
graduated to 20mm  
for hand surgery




**CASTROVIEJO CALIPER**

OM 53-0042  
6 1/2" (165mm)  
graduated to 40mm




**POLK FINGER GONIOMETER**

OM 48-0310  
6" (152mm)  
cm and inch graduated



**ALUMINUM CALIPER**

OM 48-0312  
9 1/8" (232mm)  
inch and mm graduated  
up to 6" (180mm)  
inside and outside measuring  
faces and depth bar



Inside measuring faced  
extra long (60mm) for  
neurosurgery.

**NEURO-CALIPER**

OM 48-0315  
9" (229mm)  
stainless steel  
inch and mm graduated  
up to 5 3/4" (145mm)  
inside and outside measuring  
faces and depth bar



OM 39-0692 6" (152mm)  
 OM 39-0694 8" (203mm)  
 OM 39-0698 12" (305mm)  
 OM 39-0700 20" (508mm)

**RULER**


flexible  
 inch and mm graduated



Ruler for x-rays.


**X-RAY SCALE**

OM 48-0299  
 50cm length  
 mm graduated



**MEASURING TAPE**

OM 48-0313  
 linen, inch and  
 cm graduated  
 up to 150cm



Front Back

Reverse side has handy gauge to  
 measure K-Wires, Steinmann Pins,  
 twist drills, and pins and rods.

**Note:**

Gauge Measures	
K-Wires	.028"- .062" .7mm-1.6mm
Pins	5/64"-3/16" 2.0mm-4.76mm
Rods	.197"-1/4" 5.0mm-6.35mm

**K-Wire RULER & PIN GAUGE**

OM 39-0710  
 6" (152mm)  
 inch and mm graduated



**EASY CLEAN DEPTH GAUGE**

	Screw Size	Hook Probe	Measuring Range
OM 24-0116	1.2-1.5mm	0.8mm	0-30mm
OM 24-0118	2.0-2.4mm	1.2mm	0-30mm
OM 24-0120	2.0-2.4mm	1.2mm	0-50mm
OM 24-0124	2.7mm	1.7mm	0-50mm
OM 24-0128	3.5-4.0mm	2.2mm	0-50mm
OM 24-0130	3.5-4.0mm	2.2mm	0-80mm