

CAST ROOM INSTRUMENTS



To cut or break bandages, stockinette, drapes, felt, plaster or aluminum casts.

SECTION E

CONTENTS



Scissors and Shears	1-3
Saws	3
Knives	3
Cast Breakers	3-4
Cast Spreaders	4-5

SUPER-CUT BANDAGE SCISSORS

OrthoMed developed these remarkably sharp scissors in response to the demands of doctors, nurses, and other O.R. personnel.

Hospital tests show that these scissors cut multiple layers of stockinette, bandages, drapes, and moleskin with ease.



- Very little hand strength is required for cutting, putting less strain on thumb and forearm.
- Perfect cutting control because you don't have to force them to cut.
- Scissors cut faster, saving critical time in the O.R.
- Retain their sharpness longer than ordinary scissors.
- Not recommended for cutting plaster.




All the benefits of our Super-Cut Lister scissors plus an oversized finger ring, which helps relieve pressure on the third metacarpal.

SUPER-CUT LISTER BANDAGE SCISSORS WITH LARGE FINGER RING

OM 35-0214
8" (203mm)



OM 35-0210



OM 35-0212


Cuts bandages and multiple layers of dressings with ease. Razor-sharp blades reduce hand fatigue.

SUPER-CUT LISTER BANDAGE SCISSORS

OM 35-0210
5 1/2" (140mm)
stainless

OM 35-0212
7 1/4" (184mm)
stainless


PLASTER CAST SCISSORS AND SHEARS



For light plaster and cloth.

UTILITY SCISSORS

OM 31-0832 6" (152mm) Lister-type (not pictured)
 OM 31-0830 7 1/2" (191mm)
 one serrated blade
 with autoclavable plastic finger rings



For light weight plaster or aluminum casts.

HERCULES SCISSORS


OM 31-0834
 7 1/2" (191mm)
 one serrated blade



Effective for cloth or cotton wadding.

ORTHOMED UTILITY SHEARS


OM 31-0836
 8" (203mm)
 one serrated blade
 with wire cutting notch



For heavier plaster.

BRUNS SHEARS


OM 46-0090
 9 1/4" (235mm)
 one serrated blade



Useful for cutting tape and bandages.

ESMARCH SHEARS

OM 46-0060
8" (203mm)



ENGEL PLASTER SAW

OM 46-0056
6" (152mm)
with stainless blade
chrome handle



REINER PLASTER KNIFE

OM 46-0058
7" (178mm)
with stainless handle




For lightweight cast breaking and trimming.



PEDIATRIC CAST BREAKER

OM 46-0188
7" (178mm)
deep serrated tips

PLASTER CAST INSTRUMENTS



OM 46-0186


OM 46-0184

For cast removal.

PLASTER BREAKER WOLFF (BOEHLER)



deep serrated tips
OM 46-0186
7" (178mm)
small

OM 46-0184
9 1/2" (241mm)



ORTHOMED CAST SPREADER

OM 46-0134
9" (229mm) serrated
3 prongs





*To relieve cast pressure
or for cast removal. One
handed operation.*

WALTON CAST SPREADER

OM 46-0133
9" (229mm) serrated
3 prongs


*To relieve cast pressure
or for cast removal.*



**HENNIG PLASTER
SPREADER**

OM 46-0138
11" (279mm)
serrated

For cast removal.



STILLE PLASTER SHEARS

OM 46-0156
15" (381mm)