

BONE FILES AND RASPS



To smooth, sculpt and clean bone.

SECTION F


CONTENTS

Bone Files	1-7
Bone Rasps	1-7




BONE FILE #92B

OM 69-0534
 6" (152mm)
 single ended
 2mm width angled
 plain serrations
 delicate




BONE FILE #92A

OM 69-0530
 6" (152mm)
 single ended
 2mm width angled
 plain serrations
 delicate



NAIL RASP #93


OM 69-0538
 6 3/4" (172mm)
 double ended
 2mm width angled




BONE RASP #64

OM 58-0070
 7" (178mm)
 double ended
 3mm and 4mm width

BONE FILES AND RASPS



OM 69-0541



OM 69-0540

POLOKOFF RASP

Plain serrations

OM 69-0541

7" (178mm)

double ended


3mm and 4mm width

OM 69-0540

8 1/2" (216mm)

double ended

6mm width



BONE FILE # 10



OM 69-0542

7" (178mm)

double ended

5mm width

plain serrations



BONE FILE # 12

OM 69-0543


7" (178mm)

double ended


7mm and 5mm width

plain serrations







BONE FILE # 12A
 OM 69-0544
 7" (178mm)
 double ended
 5mm width
 plain and cross serrations



BONE FILE # 12CA
 OM 69-0545
 7" (178mm)
 double ended
 5mm and 7mm width
 cross serrations




BONE FILE # 33
 OM 69-0546
 7" (178mm)
 double ended
 6mm and 5mm width
 plain serrations




BONE FILE # 45
 OM 69-0549
 7 1/4" (178mm)
 double ended
 6mm and 5mm width
 plain serrations

BONE FILES AND RASPS




For rapid shaping of bone.




BONE RASP

OM 69-0547
8 1/2" (216mm)
double ended
13mm width
fine and coarse cross serrations




JOSEPH RASP

OM 58-0068
6 1/4" (159mm)
single ended
8mm at widest point
fine cross serrations



OM 58-0072




OM 58-0073

LEWIS RASP

Plain serrations
OM 58-0072
7 1/2" (191mm)
fine
single ended
8mm width

OM 58-0073
7 1/2" (191mm)
coarse
single ended
8mm width



COTTLE RASP

OM 58-0066
8" (203mm)
single ended
backward cutting
7.5mm width
plain serrations




MALTZ RASP

OM 58-0060
7" (178mm) single ended
backward cutting
10mm width



MALTZ RASP IMPROVED

OM 58-0062
8 1/4" (210mm)
double ended
8mm width
forward and backward cutting



OM 58-0050 OM 58-0052

AUFRICHT RASP

OM 58-0050
8 1/4" (210mm)
single ended
forward cutting
9mm width


OM 58-0052
8 1/4" (210mm)
single ended
backward cutting
9mm width



FOMON RASP


OM 46-0655
8 1/2" (216mm)
double ended
fine and coarse teeth
4-sided, convex and flat blades
8mm width

BONE FILES AND RASPS




KLEINERT KUTZ RASP

OM 46-0653
 8 1/2" (216mm)
 double ended
 fine and coarse teeth
 4-sided, convex and flat blades
 8mm width




BONE FILE

OM 46-0650
 9 1/2" (241mm)
 single ended
 bayonet shank
 half round blade
 fine teeth
 15mm width




PUTTI BONE RASP

OM 46-0654
 10 1/2" (267mm)
 double ended
 flat blades
 18mm width tapers to 4mm
 one end curved up
 one end curved sideways



PUTTI BONE RASP

OM 46-0652
 12" (305mm)
 double ended
 one round blade
 10mm width tapers to 4mm
 and one half round blade
 15mm width tapers to 5mm



BONE FILE

OM 46-0705
11" (280mm)
single ended
one round blade
2-sided, convex and flat blade
23mm width tapers to 6mm