

BONE CURETTES




To scrape, shape and clean bone.

SECTION H

CONTENTS

Bone Curettes	1-6
Endoscopic Bone Curettes	7-8




JANSEN CURETTE

OM 29-0632
6" (152mm)
double ended
oval cups



MARTINI BONE CURETTE

OM 29-0620
5 1/2" (140mm)
double ended
4mm and 5mm round cups



WILLIGER BONE CURETTE


OM 29-0630
5 1/2" (140mm)
double ended
3mm and 4mm oval cups




VOLKMANN CURETTE

Oval and round cups
OM 46-0560 5" double ended
OM 46-0562 6 1/2" double ended
OM 46-0561 8" double ended

BONE CURETTES



OM 49-0200 OM 49-0202



OM 46-0550 OM 49-0204

VOLKMANN CURETTE


Oval cups, double ended

OM 49-0200 5" (127mm)


OM 49-0202 5 1/2" (140mm)

OM 46-0550 6 1/2" (165mm)

OM 49-0204 8" (203mm)



3mm
4mm
5mm



3mm
4mm
5mm



3mm
4mm
5mm



3mm
4mm
5mm

Vertical angle 45°

OM 49-0403 3mm diameter
OM 49-0404 4mm diameter
OM 49-0405 5mm diameter

Vertical angle 90°

OM 49-0413 3mm diameter
OM 49-0414 4mm diameter
OM 49-0415 5mm diameter

Horizontal angle 45°


OM 49-0423 3mm diameter
OM 49-0424 4mm diameter
OM 49-0425 5mm diameter

Horizontal angle 90°

OM 49-0433 3mm diameter
OM 49-0434 4mm diameter
OM 49-0435 5mm diameter

RAY PITUITARY CURETTES

7 1/2"




EPSTEIN CURETTE


With hollow handle
 OM 49-0580 8" (203mm) #000
 OM 49-0582 8" (203mm) #00
 OM 49-0584 8" (203mm) #0



OM 49-0574
 straight






OM 49-0576
 down angle



OM 49-0578
 up angle

SCOVILLE CURETTE

With hollow handle
 10" (254mm)





6 3/4" (172mm) straight
 OM 49-0120 #3/0
 OM 49-0121 #2/0
 OM 49-0122 #0
 OM 49-0123 #1
 OM 49-0124 #2
 OM 49-0125 #3
 OM 49-0126 #4

6 3/4" (172mm) angled 45°
 OM 49-0120A #3/0
 OM 49-0121A #2/0
 OM 49-0122A #0
 OM 49-0123A #1
 OM 49-0124A #2
 OM 49-0125A #3
 OM 49-0126A #4

BRUN (SPRATT) CURETTE

With hollow handle, round cup

BONE CURETTES



7" (178mm) straight

- OM 46-0511 #5/0
- OM 46-0512 #4/0
- OM 46-0513 #3/0
- OM 46-0515 #2/0
- OM 46-0517 #0
- OM 46-0519 #1
- OM 46-0521 #2
- OM 46-0523 #3
- OM 46-0525 #4
- OM 46-0529 #5
- OM 46-0531 #6

7" (178mm) angled 45°

- OM 46-0511A #5/0
- OM 46-0512A #4/0
- OM 46-0513A #3/0
- OM 46-0515A #2/0
- OM 46-0517A #0
- OM 46-0519A #1
- OM 46-0521A #2
- OM 46-0523A #3
- OM 46-0525A #4
- OM 46-0529A #5
- OM 46-0531A #6

9" (229mm) straight

- OM 46-0563 #3/0
- OM 46-0564 #2/0
- OM 46-0565 #0
- OM 46-0566 #1
- OM 46-0567 #2
- OM 46-0568 #3
- OM 46-0569 #4
- OM 46-0570 #5
- OM 46-0571 #6



9" (229mm) angled 45°

- OM 49-0009 #4/0
- OM 49-0110 #3/0
- OM 49-0111 #2/0
- OM 49-0112 #0
- OM 49-0113 #1
- OM 49-0114 #2
- OM 49-0115 #3
- OM 49-0116 #4
- OM 49-0117 #5
- OM 49-0118 #6

BRUN (SPRATT) CURETTE

With hollow handle, oval cup

See Section AA for Brun Rack.



Straight

- OM 46-0551 #3/0
- OM 46-0552 #2/0
- OM 46-0553 #0
- OM 46-0554 #1
- OM 46-0555 #2
- OM 46-0556 #3
- OM 46-0557 #4
- OM 46-0558 #5
- OM 46-0559 #6


Angled 45°

- OM 49-0150 #3/0
- OM 49-0151 #2/0
- OM 49-0152 #0
- OM 49-0153 #1
- OM 49-0154 #2
- OM 49-0155 #3
- OM 49-0156 #4
- OM 49-0157 #5
- OM 49-0158 #6

BRUN CURETTE

9" (229mm) with solid hexagonal handle

See Section AA for Brun Rack.



Anodized aluminum handle
stainless steel, working end

Straight

- OM 46-0422 #3/0
- OM 46-0423 #2/0
- OM 46-0424 #0
- OM 46-0425 #1
- OM 46-0426 #2
- OM 46-0427 #3
- OM 46-0428 #4
- OM 46-0429 #5
- OM 46-0430 #6

Angled

- OM 46-0433 #3/0
- OM 46-0434 #2/0
- OM 46-0435 #0
- OM 46-0436 #1
- OM 46-0437 #2
- OM 46-0438 #3
- OM 46-0439 #4
- OM 46-0440 #5
- OM 46-0441 #6

All stainless,
Straight

- OM 46-0422S #3/0
- OM 46-0423S #2/0
- OM 46-0424S #0
- OM 46-0425S #1
- OM 46-0426S #2
- OM 46-0427S #3
- OM 46-0428S #4
- OM 46-0429S #5
- OM 46-0430S #6


Angled

- OM 46-0433S #3/0
- OM 46-0434S #2/0
- OM 46-0435S #0
- OM 46-0436S #1
- OM 46-0437S #2
- OM 46-0438S #3
- OM 46-0439S #4
- OM 46-0440S #5
- OM 46-0441S #6

COBB CURETTE

11" (280mm) with knurled handle

See Section AA for Cobb Rack.



0

- OM 55-0680 10" (250mm)
- OM 55-0684 11 1/2" (292mm)



2

- OM 55-0681 10" (250mm)
- OM 55-0685 11 1/2" (292mm)




6

- OM 55-0682 10 1/4" (260mm)
- OM 55-0686 12" (305mm)

RING CURETTE

Hollow handle with fenestrated cup

BONE CURETTES



Straight


OM 49-0450 size #1, 8mm diameter
 OM 49-0452 size #2, 6mm diameter
 OM 49-0454 size #3, 3mm diameter

Angled 35°

OM 49-0456 size #1, 8mm diameter
 OM 49-0458 size #2, 6mm diameter
 OM 49-0460 size #3, 3mm diameter

CONE RING CURETTE

9" (229mm) with knurled handle
 with fenestrated cup



PIFFARD CURETTE

5 1/2" (140mm)
 OM 29-0503 small
 OM 29-0505 medium
 OM 29-0507 large


FOX CURETTE

5 1/2" (140mm) round
 OM 29-0495 1mm
 OM 29-0496 2mm
 OM 29-0497 3mm
 OM 29-0498 4mm
 OM 29-0499 5mm
 OM 29-0500 6mm


Phenolic handles are autoclavable to 290°F (143°C).




ENDOSCOPIC BRUN CURETTE
Graduated with phenolic handle




OM 55-0532 17" (432mm) 40° angle
2mm x 4mm cup
2.7" tip length
9.3" shaft length
5" handle length
.285" (7.2mm) shaft diameter




OM 55-0523 17" (432mm) straight
3mm x 6mm cup
2.7" tip length
9.3" shaft length
5" handle length
.285" (7.2mm) shaft diameter



OM 55-0533 17" (432mm) 40° angle
3mm x 6mm cup
2.7" tip length
9.3" shaft length
5" handle length
.285" (7.2mm) shaft diameter




OM 55-0543 20" (508mm) 40° angle
3mm x 6mm cup
3.5" tip length
10" shaft length
6.5" handle length
.45" (11.5mm) shaft diameter




OM 55-0556 20" (508mm) 30° angle
6mm x 9mm cup
3.5" tip length
10" shaft length
6.5" handle length
.45" (11.5mm) shaft diameter



ENDOSCOPIC EPSTEIN CURETTE
Graduated with phenolic handle





OM 55-0593 17" (432mm)
reverse 40° angle
3mm x 6mm cup
2.7" tip length
9.3" shaft length
5" handle length
.285" (7.2mm) shaft diameter



OM 55-0586 20" (508mm)
reverse 30° angle
6mm x 9mm cup
3.5" tip length
10" shaft length
6.5" handle length
.45" (11.5mm) shaft diameter

Phenolic handles are autoclavable to 290°F (143°C).


ENDOSCOPIC BONE CURETTES




OM 55-0510 20" (508mm)
 straight
 9mm x 15mm cup
 3.5" tip length
 10" shaft length
 6.5" handle length
 .45" (11.5mm) shaft diameter

Phenolic handles are autoclavable to 290°F (143°C).

ENDOSCOPIC SCOVILLE CURETTE
 With phenolic handle



OM 55-0573 17" (432mm)
 40° angle
 3mm cup diameter
 2.7" tip length
 9.3" shaft length
 5" handle length
 .285" (7.2mm) shaft diameter



OM 55-0566 20" (508mm)
 30° angle
 6mm cup diameter
 3.5" tip length
 10" shaft length
 6.5" handle length
 .45" (11.5mm) shaft diameter

Phenolic handles are autoclavable to 290°F (143°C).

ENDOSCOPIC CONE RING CURETTE
 With phenolic handle



LONG BONE CURETTE

Straight 15" (381mm)

- OM PS00 991-01 #4/0
- OM PS00 991-02 #3/0
- OM PS00 991-03 #2/0
- OM PS00 991-04 #0
- OM PS00 991-05 #1
- OM PS00 991-06 #2
- OM PS00 991-07 #3
- OM PS00 991-08 #4
- OM PS00 991-09 #5
- OM PS00 991-10 #6

LONG BONE CURETTE

Angled 15" (381mm)

- OM PS00 991-01A #4/0
- OM PS00 991-02A #3/0
- OM PS00 991-03A #2/0
- OM PS00 991-04A #0
- OM PS00 991-05A #1
- OM PS00 991-06A #2
- OM PS00 991-07A #3
- OM PS00 991-08A #4
- OM PS00 991-09A #5
- OM PS00 991-10A #6



VOLKMANN CURETTE

Straight 15" (381mm)

- OM 55-0576 (10mm X 19mm cup)
- OM 55-0577 (13mm X 25mm cup)
- OM 55-0578 (15mm X 29mm cup)



LONG COBB CURETTE

15" overall length with 9" handle

Straight

OM 93-0112	15" (381mm) overall, #3/0
OM 93-0113	15" (381mm) overall, #2/0
OM 93-0114	15" (381mm) overall, #0
OM 93-0115	15" (381mm) overall, #1
OM 93-0116	15" (381mm) overall, #2
OM 93-0117	15" (381mm) overall, #3
OM 93-0118	15" (381mm) overall, #4
OM 93-0119	15" (381mm) overall, #5
OM 93-0125	15" (381mm) overall, #6
OM 93-0127	15" (381mm) overall, #7

Angled

OM 93-0205	15" (381mm) overall, #3/0
OM 93-0206	15" (381mm) overall, #2/0
OM 93-0207	15" (381mm) overall, #0
OM 93-0208	15" (381mm) overall, #1
OM 93-0209	15" (381mm) overall, #2
OM 93-0210	15" (381mm) overall, #3
OM 93-0211	15" (381mm) overall, #4
OM 93-0212	15" (381mm) overall, #5
OM 93-0213	15" (381mm) overall, #6
OM 93-0214	15" (381mm) overall, #7