


OSTEOTOMES, CHISELS AND GOUGES



To score, cut, scrape, clean and sculpt bone.

SECTION J


CONTENTS

Interchangeable Set	1
Chisels	4-9
Gouges	5-2
Osteotomes	2-12



OSTEOTOME BLADES

- OM 49-0668 5mm
- OM 49-0672 10 mm
- OM 49-0676 16 mm
- OM 49-0680 25 mm




OM 49-0672




CHISEL BLADES

- OM 49-0652 5mm
- OM 49-0656 10 mm
- OM 49-0660 16 mm
- OM 49-0664 25 mm




OM 49-0664



KEY

OM 49-0634



OM 49-0640

OM 49-0644

OM 49-0648

GOUGE BLADES

- OM 49-0640 1" width, 50mm radius
- OM 49-0644 1" width, 60 mm radius
- OM 49-0648 1" width, 70 mm radius

Blades should be discarded after use.

Phenolic handle is autoclavable to 290°F (143°C).

INTERCHANGEABLE OSTEOTOME, CHISEL AND GOUGE SET

- OM 49-0630 complete set with handle, key and all blades
- OM 49-0632 7" phenolic handpiece only

OSTEOTOMES, CHISELS AND GOUGES

5" Straight

OM 46-0336	4mm	OM 46-0343	14mm
OM 46-0338	6mm	OM 46-0345	15mm
OM 46-0340	8mm	OM 46-0346	18mm
OM 46-0342	10mm	OM 46-0347	19mm
OM 46-0344	12mm	OM 46-0335	20mm

MINI-LAMBOTTE OSTEOTOME

7" Straight
With calibration lines

OM 69-0940	10mm
OM 69-0946	20mm
OM 69-0950	25mm

7" Curved
With calibration lines

OM 69-0970	5mm
OM 69-0975	10mm

The unusually long bevel allows the cut bone to bend gradually as it is cut.

LONG BEVEL OSTEOTOME


7" Straight
With calibration lines

OM 46-0348	4mm
OM 46-0349	5mm
OM 46-0350	6mm
OM 46-0358	8mm
OM 46-0351	10mm
OM 46-0352	12mm
OM 46-0353	16mm
OM 46-0354	19mm
OM 46-0355	22mm
OM 46-0356	25mm

7" Curved
With calibration lines

OM 46-0348C	4mm
OM 46-0349C	5mm
OM 46-0350C	6mm
OM 46-0358C	8mm
OM 46-0351C	10mm
OM 46-0352C	12mm
OM 46-0353C	16mm
OM 46-0354C	19mm
OM 46-0355C	22mm
OM 46-0356C	25mm

LAMBOTTE OSTEOTOME



9" Straight
Without calibration lines


OM 46-0304	1/4"	6mm
OM 46-0306	1/2"	13mm
OM 46-0310	3/4"	19mm
OM 46-0314	1"	25mm
OM 46-0318	1 1/4"	32mm
OM 46-0322	1 1/2"	38mm

9" Curved
Without calibration lines

OM 46-0328	1/4"	6mm
OM 46-0329	1/2"	13mm
OM 46-0330	3/4"	19mm
OM 46-0331	1"	25mm
OM 46-0332	1 1/4"	32mm
OM 46-0333	1 1/2"	38mm

LAMBOTTE OSTEOTOME

Racks for 9" Lambotte Osteotomes are shown on page AA2.



Straight

OM 46-0497	1/8"	3mm
OM 46-0498	3/16"	4mm
OM 46-0499	1/4"	6mm
OM 46-0500	5/16"	8mm
OM 46-0501	3/8"	10mm
OM 46-0502	1/2"	13mm
OM 46-0503	5/8"	16mm
OM 46-0504	3/4"	19mm



Curved

OM 46-0700	1/8"	3mm
OM 46-0701	3/16"	4mm
OM 46-0702	1/4"	6mm
OM 46-0703	5/16"	8mm
OM 46-0704	3/8"	10mm
OM 46-0706	1/2"	13mm
OM 46-0708	5/8"	16mm
OM 46-0709	3/4"	19mm

HOKE OSTEOTOME
5 1/4" (140mm) straight

Racks for Hoke Osteotomes are shown on page AA3.

OSTEOTOMES, CHISELS AND GOUGES



- OM 46-0505 3/16" (4mm)
- OM 46-0506 1/4" (6mm)
- OM 46-0507 3/8" (10mm)
- OM 46-0508 1/2" (13mm)
- OM 46-0509 3/4" (19mm)

HOKE CHISEL

5 1/4" (133mm)
straight


Racks for Hoke Chisels
are shown in Section AA.




- OM 46-0268 2mm
- OM 46-0270 4mm
- OM 46-0272 6mm
- OM 46-0274 8mm
- OM 46-0276 10mm
- OM 46-0278 12mm
- OM 46-0280 15mm

SHEEHAN OSTEOTOME

6 1/4" (159mm)
straight



OM 58-0102





OM 58-0104

- OM 58-0100 straight
- OM 58-0102 curved, guard right
- OM 58-0104 curved, guard left

SILVER OSTEOTOME




7" (178mm)
with single guard




- OM 46-0284 6mm
- OM 46-0286 12mm
- OM 46-0288 18mm

U.S. ARMY PATTERN OSTEOTOME

7" (178mm)



6mm



12mm

U.S. ARMY PATTERN CHISEL
7" (178mm)


OM 46-0289	6mm
OM 46-0290	12mm
OM 46-0291	14mm

6mm






12mm




14mm

U.S. ARMY PATTERN GOUGE
7" (178mm)


OM 46-0292	6mm
OM 46-0293	12mm
OM 46-0294	14mm


4mm




6mm




8mm



10mm






12mm




14mm

ALEXANDER CHISEL
7" (178mm)


OM 49-0040	4mm
OM 49-0041	6mm
OM 49-0042	8mm
OM 49-0043	10mm
OM 49-0044	12mm
OM 49-0045	14mm


4mm




6mm




8mm



10mm



12mm




14mm

ALEXANDER GOUGE
7" (178mm)


OM 49-0050	4mm
OM 49-0051	6mm
OM 49-0052	8mm
OM 49-0053	10mm
OM 49-0054	12mm
OM 49-0055	14mm

OSTEOTOMES, CHISELS AND GOUGES



OM 49-0260	4mm
OM 49-0262	6mm
OM 49-0264	8mm
OM 49-0266	10mm
OM 49-0268	12mm

LEXER MINI CHISEL
7" (178mm)



OM 49-0246	5mm
OM 49-0248	7mm
OM 49-0250	10mm
OM 49-0252	15mm
OM 49-0254	20mm
OM 49-0256	25mm
OM 49-0258	30mm

LEXER CHISEL
8 3/4" (222mm)

OSTEOTOMES, CHISELS AND GOUGES



1/4"



1/2"



1"

8" Straight

OM 46-0448	1/4"	6mm
OM 46-0449	3/8"	10mm
OM 46-0450	1/2"	13mm
OM 46-0451	5/8"	16mm
OM 46-0452	3/4"	19mm
OM 46-0453	1"	25mm
OM 46-0454	1 1/4"	32mm


8" Curved

OM 46-0457	1/4"	6mm
OM 46-0458	3/8"	10mm
OM 46-0459	1/2"	13mm
OM 46-0460	5/8"	16mm
OM 46-0461	3/4"	19mm
OM 46-0462	1"	25mm
OM 46-0463	1 1/4"	32mm




SMITH PETERSON OSTEOTOME


Racks for Smith Peterson Osteotomes are shown in Section AA.



1/4"



1/2"



3/4"

8" Straight

OM 46-0466	1/4"	6mm
OM 46-0467	3/8"	10mm
OM 46-0468	1/2"	13mm
OM 46-0476	5/8"	16mm
OM 46-0477	3/4"	19mm
OM 46-0478	1"	25mm
OM 46-0479	1 1/4"	32mm

8" Curved

OM 46-0487	1/4"	6mm
OM 46-0488	3/8"	10mm
OM 46-0489	1/2"	13mm
OM 46-0490	5/8"	16mm
OM 46-0491	3/4"	19mm
OM 46-0493	1"	25mm
OM 46-0495	1 1/4"	32mm



SMITH PETERSON GOUGE

Racks for Smith Peterson Gouges are shown in Section AA.

OSTEOTOMES, CHISELS AND GOUGES



- OM 49-0002 10mm
- OM 49-0004 14mm
- OM 49-0006 18mm
- OM 49-0008 22mm

STILLE-TYPE OSTEOTOME
8" (203mm) straight





- OM 49-0010 10mm
- OM 49-0012 14mm
- OM 49-0014 18mm
- OM 49-0016 22mm

STILLE-TYPE CHISEL
8" (203mm) straight



- OM 49-0018 10mm
- OM 49-0020 14mm
- OM 49-0022 18mm
- OM 49-0024 22mm

STILLE-TYPE GOUGE
8" (203mm) straight



9" (229mm) Straight


- OM 46-0359 1/4" (6mm)
- OM 46-0360 3/8" (10mm)
- OM 46-0361 1/2" (13mm)
- OM 46-0362 5/8" (16mm)
- OM 46-0363 3/4" (19mm)
- OM 46-0364 7/8" (22mm)
- OM 46-0365 1" (25mm)
- OM 46-0366 1 1/8" (29mm)
- OM 46-0367 1 1/4" (31mm)
- OM 46-0368 1 1/2" (38mm)

9" (229mm) Curved

- OM 46-0373 1/4" (6mm)
- OM 46-0376 3/8" (10mm)
- OM 46-0377 1/2" (13mm)
- OM 46-0378 5/8" (16mm)
- OM 46-0379 3/4" (19mm)
- OM 46-0380 7/8" (22mm)
- OM 46-0381 1" (25mm)
- OM 46-0382 1 1/8" (29mm)
- OM 46-0383 1 1/4" (31mm)
- OM 46-0384 1 1/2" (38mm)

HIBBS OSTEOTOME

Racks for Hibbs Osteotomes are shown in Section AA.




1/4"




1/2"



3/4"



1"



1 1/4"

9" (229mm) Straight

OM 49-0059	1/4" (6mm)
OM 49-0060	3/8" (10mm)
OM 49-0061	1/2" (13mm)
OM 49-0062	5/8" (16mm)
OM 49-0063	3/4" (19mm)
OM 49-0064	7/8" (22mm)
OM 49-0065	1" (25mm)
OM 49-0066	1 1/8" (29mm)
OM 49-0067	1 1/4" (31mm)
OM 49-0068	1 1/2" (38mm)

9" (229mm) Curved

OM 49-0073	1/4" (6mm)
OM 49-0074	3/8" (10mm)
OM 49-0075	1/2" (13mm)
OM 49-0076	5/8" (16mm)
OM 49-0077	3/4" (19mm)
OM 49-0078	7/8" (22mm)
OM 49-0079	1" (25mm)
OM 49-0080	1 1/8" (29mm)
OM 49-0081	1 1/4" (31mm)
OM 49-0082	1 1/2" (38mm)

HIBBS CHISEL


Racks for Hibbs Chisel are shown in Section AA.




1/4"



1/2"



3/4"



1"



1 1/4"



9" (229mm) Straight

OM 46-0387	1/4" (6mm)
OM 46-0388	3/8" (10mm)
OM 46-0389	1/2" (13mm)
OM 46-0390	5/8" (16mm)
OM 46-0391	3/4" (19mm)
OM 46-0392	7/8" (22mm)
OM 46-0393	1" (25mm)
OM 46-0394	1 1/8" (29mm)
OM 46-0395	1 1/4" (31mm)
OM 46-0401	1 1/2" (38mm)



9" (229mm) Curved

OM 46-0402	1/4" (6mm)
OM 46-0403	3/8" (10mm)
OM 46-0404	1/2" (13mm)
OM 46-0405	5/8" (16mm)
OM 46-0407	3/4" (19mm)
OM 46-0408	7/8" (22mm)
OM 46-0409	1" (25mm)
OM 46-0410	1 1/8" (29mm)
OM 46-0411	1 1/4" (31mm)
OM 46-0412	1 1/2" (38mm)

HIBBS GOUGE



Racks for Hibbs Gouge are shown in Section AA.

OSTEOTOMES, CHISELS AND GOUGES






3mm

SHOULDER PENETRATING GOUGE
 OM 46-0919
 8 1/2" (220mm)

SWAN NECK GOUGE
 OM 49-0084
 9" (229mm)

11" (280mm) Straight		11" (280mm) Curved	
OM 49-0350	1/4" (6mm)	OM 49-0357	1/4" (6mm)
OM 49-0351	3/8" (10mm)	OM 49-0358	3/8" (10mm)
OM 49-0352	1/2" (13mm)	OM 49-0359	1/2" (13mm)
OM 49-0353	5/8" (16mm)	OM 49-0360	5/8" (16mm)
OM 49-0354	3/4" (19mm)	OM 49-0361	3/4" (19mm)
OM 49-0355	1" (25mm)	OM 49-0362	1" (25mm)
OM 49-0356	1 1/4" (31mm)	OM 49-0363	1 1/4" (31mm)


COBB OSTEOTOME



OM 49-0030
straight
5/32" (7mm) blade width



OM 49-0032
medium curve
5/32" (7mm) blade width




OM 49-0034
strong curve
5/32" (7mm) blade width



OM 49-0036
reverse curve
5/32" (7mm) blade width


COBB-TYPE SPINAL GOUGE
11" (280mm)
All stainless, with knurled handle

Racks for Cobb-Type Gouges
are shown in Section AA.



ENDOSCOPIC HIBBS GOUGE

OM 55-0630
18 1/2" (470mm) with stainless steel handle
.28" (7mm) tip width straight
3 1/2" tip length
10 3/4" shaft length
4 1/4" handle length
.45" (11.5mm) shaft diameter



ENDOSCOPIC HIBBS OSTEOTOME

OM 55-0642
18 1/2" (470mm) with stainless steel handle
3/8" (10mm) tip width straight
3 1/2" tip length
10 3/4" shaft length
4 1/4" handle length
.45" (11.5mm) shaft diameter