

MEASURING DEVICES



To determine dimensions during surgery.

SECTION W


CONTENTS

Screw Depth Gauges	1
Calipers	1-2
Rulers	3



TOWNLEY CALIPER

OM 48-0314
4" (102mm)
takes inside and outside measurements up to 4" (10cm)



BONE COMPASS

OM 48-0302
8" (203mm)

Used to determine screw length required during internal fracture fixation (mm graduated). Gauges snap apart for easy cleaning.



MINI SCREW DEPTH GAUGE

OM 47-0934
4 3/4" (121mm)
up to 30mm



SMALL SCREW DEPTH GAUGE

OM 47-0936
6 1/4" (159mm)
up to 50mm




LARGE SCREW DEPTH GAUGE

OM 47-0940
8 3/4" (222mm)
up to 100mm




X-LARGE SCREW DEPTH GAUGE

OM 47-0935
11" (279mm)
up to 120mm




CASTROVIEJO CALIPER

OM 53-0040
 3 1/2" (89mm)
 graduated to 20mm
 for hand surgery




CASTROVIEJO CALIPER

OM 53-0042
 6 1/2" (165mm)
 graduated to 40mm




POLK FINGER GONIOMETER

OM 48-0310
 6" (152mm)
 cm and inch graduated



ALUMINUM CALIPER

OM 48-0312
 9 1/8" (232mm)
 inch and mm graduated
 up to 6" (180mm)
 inside and outside measuring
 faces and depth bar



*Inside measuring faced
 extra long (60mm) for
 neurosurgery.*


NEURO-CALIPER

OM 48-0315
 9" (229mm)
 stainless steel
 inch and mm graduated
 up to 5 3/4" (145mm)
 inside and outside measuring
 faces and depth bar




OM 39-0692 6" (152mm)
 OM 39-0694 8" (203mm)
 OM 39-0698 12" (305mm)
 OM 39-0700 20" (508mm)

RULER
 flexible
 inch and mm graduated




Ruler for x-rays.

X-RAY SCALE
 OM 48-0299
 50cm length
 mm graduated



MEASURING TAPE
 OM 48-0313
 linen, inch and
 cm graduated
 up to 150cm



Front Back

*Reverse side has handy gauge to
 measure K-Wires, Steinmann Pins,
 twist drills, and pins and rods.*

Note:
 Gauge Measures
 K-Wires .028"-.062"
 .7mm-1.6mm
 Pins 5/64"-3/16"
 2.0mm-4.76mm
 Rods .197"-1/4"
 5.0mm-6.35mm

K-WIRE RULER & PIN GAUGE
 OM 39-0710
 6" (152mm)
 inch and mm graduated