


MISCELLANEOUS INSTRUMENTS



SECTION Z

CONTENTS

Lead Hands	1
Forceps Holders	2
Finger Ring Cutter	2
Pin Wheel	2

Note: Oxidation of lead surfaces is a normal process.



LEAD HAND

OM 49-0092
9 1/2" (241mm)
child



LEAD HAND

OM 49-0090
14" (356mm)
adult


To prolong the life of lead hands, autoclave temperatures should not exceed 248°F (120°C). Avoid using chemical sterilization as they can cause severe corrosion.

Thinner and more flexible model. Longer fingers wrap on ends eliminating side tabs.



LEAD HAND

OM 49-0093
9 1/2" (241mm)
adult
2mm thick



ORTHOMED FORCEPS HOLDER

OM 39-0672 OM 39-0674
6" (152mm) 12" (305mm)




BUNT FORCEPS HOLDER

OM 39-0670



FINGER RING CUTTER

OM 29-0430 OM 29-0434
6 1/2" (165mm) chrome saw blade only
with saw blade



**WARTENBERG
NEUROLOGICAL PIN WHEEL**

OM 28-0514 7" (178mm)